


ARTHROJET

PULSED LAVAGE



ARTHROJET ergonomics allows total device control in all working conditions.

Standard Tip and Long Tip
with cone protection
Recommended for:
Knee arthroplastic surgery
General wound and ulcers cleaning

ARTHROJET combines efficient pulsed irrigation and a fluid-removing device in a single instrument.



ARTHROJET

PULSED LAVAGE

Tips and Accessories

Ref. n° 8000018	ARTHRO-JET Pulsed Lavage Standard Tip for knee with cone protection
Ref. n° 8000019	ARTHRO-JET Pulsed Lavage Standard Tip with cone protection
Ref. n° 8001048	ARTHRO-JET "Tibial Canal Tip"
Ref. n° 8001058	ARTHRO-JET "Femoral Canal Tip"
Ref. n° 8001068	ARTHRO-JET "Femoral Canal Brush"
Ref. n° 8001078	ARTHRO-JET "Standard Tip" for knee with cone protection
Ref. n° 8001088	ARTHRO-JET "Long Tip" for knee with cone protection
Ref. n° 8001098	ARTHRO-JET "Round Splash Shield" 30 cm
Ref. n° 8001108	ARTHRO-JET "Rectangular Splash Shield" 12,5 cm x 20 cm

Technical Information

Certification	CE 0476, UNI ISO 9001:2000, EN ISO 13485:2004, UNI CEI EN ISO 14971:2004 UNI-CEI-EN 60601-1 and enclosures, UNI CEI-EN 60601-1-2 and enclosures UNI EN 1041:2000, UNI EN ISO 980:2004, UNI EN ISO 10993-1:2004, UNI EN ISO 10993-7:1997, UNI EN ISO 10993-8:2002, EN 556:1996, EN 550:1996, UNI EN ISO 14937:2002, EN 550 and ISO 1135 (1994)
Sterilization	Process ETO EN 550 ISO 11135
Power	8 Alkaline Batteries AA
External Temperature	-5°C +60°C
Stock time	36 months
Weight	Gun with tube +/- 500 gr / batteries +/- 250 gr.
Tips available	Tips with air-blow system and cone protection Tips large surface Tips for tibial canal Tips for femoral canal Other tips on request
Release	Pistol: 12 pcs dispenser Tips: 10 pcs box
Product Patent	CE 0476



ARTHROJET
PULSED LAVAGE

